

Technical Data Sheet (valid from batch #80000)

## **KUNOS Polish No. 760**

**new** quality

**Application range** In the interior, for the post-treatment of solid wood furniture varnished in white with LIVOS

products (e.g. "bianco" optics). Do not use for cabinet interior surfaces.

**Technical qualities** Refreshes oiled surfaces and complements the color. Free from organic chemical compounds

according to DIN EN 71 part 9.

**Full declaration** Linseed oil - stand oil - natural resin ester, linseed oil - stand oil, linseed oil - wood

oil – stand oil, natural resin glycerol ester, titanium dioxide, castor oil – stand oil, silicic acid, micronized wax, alumina, orange oil, isoaliphates, dehydrated amino sugar, soya lecithin, and

drying agents free of lead (Ca, Mn, Zr).

Color White.

**Application method** Sand surface slightly (240 grain) and apply KUNOS Polish very thin with a brush. Polish in after

about 10 - 15 min. and allow to dry. If pretreated with products from other manufacturers: Test

polish in an inconspicuous spot in order to check the compatibility.

**Coverage** 1 liter (1.06 quart) undiluted is sufficient for approximately 14 m<sup>2</sup> (151 ft<sup>2</sup>), i.e. 71 ml/m<sup>2</sup>. Could

be considerably less, depending on the condition and absorbency of the surface.

**Drying time** At 23°C (73°F) and 50% relative humidity: approximately 2 to 8 hours.

**Specific gravity** Approximately 0.88 g/ml

Safety advice Ensure good ventilation during and after application. Always be sure to store and then dispose

of cleaning cloths, polishing pads, sponges, etc. that are still wet or even slightly moist with KUNOS Polish No. 760 in an airtight metal container or in water, as there is a danger of spontaneous combustion caused by the plant oil content. The liquid and applied product itself

will not combust spontaneously.

P102 Keep out of reach of children.

P301+P310 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. Immediately call a

doctor.

EUH208 Contains limonene (orange oil. May produce an allergic reaction.

Hotline for allergy sufferers: +49 (0)5825 / 8871.

**Disposal** According to locally regulated laws.

**Container sizes** 0.05 l; 0.75 l; 2.5 l; 10 l.

**Storage** Cool and dry. Shelf life: unopened at least 4 years. Close lid tightly after use.

These directions are the results of long years of research and practical testing. They are backed by our most current expert information. New findings may invalidate this information. The latest version is available on our website at www.livos.de. This data sheet is meant to serve as information and instruction. No legal liability should be interpreted from it. In case of doubt, please contact us.

March 2019